

invent-uk

EN 14175

Flow Containment & Instrumentation

CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH EN 14175, PART 3

CERTIFICATE & REPORT NO: INV/EN14175/328

DATE: 24th August 2005



Fume Cupboard Manufacturer:

Esco Micro Pte Ltd
21 Changi South Street 1
Singapore 486 777

Fume Cupboard Type:

1.34m wide ADC-4A1 bench-type
recirculatory filtration

Sash Opening Width: 1270mm

Sash Opening Height: 200mm

Average Face Velocity: 0.54m/s

Extract volume flow rate: 0.145m³/s

This is to certify that the fume cupboard described above has been type tested in accordance with EN 14175, Part 3 and resulted in performance characteristics given in the corresponding test report.

Tested and Certified by: Dr A F Bicen

Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR

Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

Handwritten signature of A. F. Bicen.