

invent-uk

NF X 15-203

Flow Containment & Instrumentation

CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH NF X 15 203

CERTIFICATE & REPORT NO: INV/BS7258/326
DATE: 15th August 2005



Enclosure Manufacturer:
Esco Micro Pte Ltd
21 Changi South Street 1
Singapore 486 777

Enclosure Type:
1.34m wide ADC-4A1 bench-type
with recirculating air filtration
Sash Opening Width: 1270mm
Sash Opening Height: 200mm
Average Face Velocity: 0.53m/s

This is to certify that the enclosure described above has been type tested in accordance with NF X 15 203 and resulted in performance characteristics given in the corresponding test report.

Certified by: Dr A F Bicen
Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR
Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

A handwritten signature in blue ink that reads 'A.F. Bicen'.